




neodisher® FA



Alkaline detergent for the automated reprocessing of thermostable and thermolabile instruments

Liquid concentrate

Fields of application:

- Automated cleaning of surgical instruments, anaesthetic equipment, surgical clogs and other medical utensils as well as baby bottles in hospitals and medical practices
- Automated cleaning of laboratory glassware in medical laboratories

Performance spectrum:

- Versatile detergent
- Good cleaning performance; reliably removes blood, protein, culture media and gelatine
- Suitable for surgical instruments, anaesthetic equipment, laboratory glassware or plastics
- Anodised aluminium must be tested first for suitability
- With titanium and titanium alloys colour changes might occur due to a change in the thickness of the colour-producing titanium oxide layer.
- Usable with any level of water hardness

Special properties:

- Excellent material compatibility
- Free of oxidising agents and surfactants

Application and dosage:

neodisher FA is used in washer disinfectors for instruments. The dosage depends among other things on the degree of soiling of the items to be washed and is 2 - 4 ml/l.

Use suitable dosing devices for dosing.

Application examples:

Automated reprocessing of surgical instruments and anaesthetic equipment with thermal disinfection:	
Cleaning with neodisher FA	2 - 4 ml/l
Neutralisation with neodisher Z oder neodisher N	1 - 2 ml/l
Final rinse with neodisher MediKlar	0,3 – 1 ml/l
Removal of blood, serum, culture media and analyses residues in medical and analytical laboratories:	
Cleaning with neodisher FA	2 - 4 ml/l
Neutralisation with neodisher Z oder neodisher N	1 - 2 ml/l
Cleaning of baby bottles in milk kitchens:	
Cleaning with neodisher FA	2 - 4 ml/l
Neutralisation with neodisher Z	1 - 2 ml/l

For avoiding water stains the use of deionised water in the final rinse is recommended. At the same time this protects anodised aluminium.

Notes on application:

- The neodisher FA-solution has to be rinsed off completely with water (preferably with deionised water).
- For professional use only
- Do not mix with other products
- Rinse out dosing system including suction hose with water before changing product
- Reprocessing should comply with all ordinances pursuant to the regulations on medical devices and should be performed with appropriate validated processes



- Please observe the reprocessing recommendations of the medical device manufacturers according to the requirements of the DIN EN ISO 17664
- The instructions given by the manufacturer of the washer disinfectant are to be observed

Technical data:

pH range	11.6 - 11.8 (2 - 4 ml/l, determined in deionised water, 20 °C)
Viscosity	< 50 mPa s (concentrate, 20 °C)
Density	approx. 1.5 g/cm ³ (20 °C)
Titration factor	0.65 (in accordance with the neodisher titration method)

Ingredients:


Ingredients according to Regulation (EC) No 648/2004 on detergents:
15 - 30 % phosphates

CE-mark:

neodisher FA complies with European guidelines for medical devices.

If a serious incident occurs with the product, report it to the manufacturer and the relevant national authority.

Storage information:

Sensitive to frost below - 15 °C. Always store at a temperature between -15 °C and 30 °C. Usable for 3 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hourglass symbol .

Hazard and precautionary statements:

For safety information see Safety Data Sheets. These are available e.g. at www.drweigert.com under the category "Service/Downloads".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Safety Data Sheet.

MB 4101/3-1

Date of issue: 11/2020