

FOB test - occult blood in the stool

Rapid immunochromatographic test for the determination of human occult blood in faeces.

The test uses the "sandwich" method with a double antibody capable of selectively identifying Faecal Occult Blood at a concentration equal to or greater than 50 ng / ml or 6 µg / g of faeces.

The accuracy of the test is not affected by the patient's diet.

It provides accurate results in minutes, facilitating timely patient treatment choices.


Simple procedure, easy to use and understand.

It does not require additional equipment.

Pack of 25 tests consisting of tube, test cassette, reagent.

For professional use

Manual in: GB, IT, FR, DE, ES PT.



Code	Description
71024531	FOB test - occult blood in the stool