

Micropipette mod. ACURA 825 from 10 - 100 µl

The Acura® line was created with the user in mind.

These tools are the product of the union of high-tech materials and extraordinary ergonomic features. In research applications, but also in routine administration, as well as offering exceptional performance, these products make pipetting more convenient and safer.

- Extraordinarily ergonomic
- Light and delicate
- Unmatched stability and metrological performance
- Adjustable tip ejector *, compatible with most tips
- Quick swift-set calibration system *: Allows for easy calibration in the laboratory
- Resistance to shocks, UV rays and autoclave washing
- CE IVD 98/79 EEC certification
- Three-year guarantee

* Socorex patent

Equipped with a precision digital display, the pipettes of this line are continuously adjustable to a very accurate volume definition.

Adjustment can be done easily with one hand while the other hand is engaged.

Snap-locks and free-spinning button removable cap ensure stability of adjustments during pipetting. Resistance to shocks, heat, chemicals and UV light is guaranteed by selected materials and careful construction.

As a result, almost no maintenance is required.

Deep cleaning and sterilization are achieved by autoclaving the fully assembled instrument (121 ° C, 20 minutes).

Each micropipette bears its own serial number and undergoes a rigorous functional check, attested by an individual QA certificate.

The 825 line includes micropipettes from 0.1 to 1000 µl.

Also available are packs (Triopack) of micropipettes in multiple packages to cover multiple ranges of volumes.



Triopack Acura® manual 825 e 835

Conf.	Contenuto (modelli micro-/macropipette)			No. cat.
C	2 µl (0.1 - 2)	10 µl (0.5 - 10)	50 µl (5 - 50)	825.TRIO.C
J	2 µl (0.1 - 2)	20 µl (2 - 20)	200 µl (20 - 200)	825.TRIO.J
M	10Y µl (1 - 10)	50 µl (5 - 50)	200 µl (20 - 200)	825.TRIO.M
O	10 µl (0.5 - 10)	100 µl (10 - 100)	1000 µl (100 - 1000)	825.TRIO.O
P	10Y µl (1 - 10)	100 µl (10 - 100)	1000 µl (100 - 1000)	825.TRIO.P
S	20 µl (2 - 20)	200 µl (20 - 200)	1000 µl (100 - 1000)	825.TRIO.S
T	50 µl (5 - 50)	200 µl (20 - 200)	1000 µl (100 - 1000)	825.TRIO.T
V	100 µl (10 - 100)	1000 µl (100 - 1000)	5 ml (0.5 - 5)	825.TRIO.V
Y	1000 µl (100 - 1000)	5 ml (0.5 - 5)	10 ml (1 - 10)	835.TRIO.Y



Code	Description	Volume ul	Division ul	Tip type	Weight gr	% Error - vol .min./medio /max	% Imprecision - vol.min./me dium/max
21825100	Micropipette mod. ACURA 825 from 10 - 100 µl	10-100	0.1	200 ul	92	<±1,5/<±1,2/<± 0,8	<1,0/<0,6/<0,2