

Micropipette mod. ACURA 825 from 0.1 - 2 μ l

The Acura® line was created with the user in mind.

These tools are the product of the union of high-tech materials and extraordinary ergonomic features. In research applications, but also in routine administration, as well as offering exceptional performance, these products make pipetting more convenient and safer.

- Extraordinarily ergonomic
- Light and delicate
- Unmatched stability and metrological performance
- Adjustable tip ejector *, compatible with most tips
- Quick swift-set calibration system *: Allows for easy calibration in the laboratory
- Resistance to shocks, UV rays and autoclave washing
- CE IVD 98/79 EEC certification
- Three-year guarantee

* Socorex patent

Equipped with a precision digital display, the pipettes of this line are continuously adjustable to a very accurate volume definition.

Adjustment can be done easily with one hand while the other hand is engaged.

Snap-locks and free-spinning button removable cap ensure stability of adjustments during pipetting. Resistance to shocks, heat, chemicals and UV light is guaranteed by selected materials and careful construction.

As a result, almost no maintenance is required.

Deep cleaning and sterilization are achieved by autoclaving the fully assembled instrument (121 ° C, 20 minutes).

Each micropipette bears its own serial number and undergoes a rigorous functional check, attested by an individual QA certificate.

The 825 line includes micropipettes from 0.1 to 1000 μ l.

Also available are packs (Triopack) of micropipettes in multiple packages to cover multiple ranges of volumes.



Triopack Acura® *manual* 825 e 835

Conf.	Contenuto (modelli micro-/macropipette)			No. cat.
C	2 µl (0.1 – 2)	10 µl (0.5 – 10)	50 µl (5 – 50)	825.TRIO.C
J	2 µl (0.1 – 2)	20 µl (2 – 20)	200 µl (20 – 200)	825.TRIO.J
M	10Y µl (1 – 10)	50 µl (5 – 50)	200 µl (20 – 200)	825.TRIO.M
O	10 µl (0.5 – 10)	100 µl (10 – 100)	1000 µl (100 – 1000)	825.TRIO.O
P	10Y µl (1 – 10)	100 µl (10 – 100)	1000 µl (100 – 1000)	825.TRIO.P
S	20 µl (2 – 20)	200 µl (20 – 200)	1000 µl (100 – 1000)	825.TRIO.S
T	50 µl (5 – 50)	200 µl (20 – 200)	1000 µl (100 – 1000)	825.TRIO.T
V	100 µl (10 – 100)	1000 µl (100 – 1000)	5 ml (0.5 – 5)	825.TRIO.V
Y	1000 µl (100 – 1000)	5 ml (0.5 – 5)	10 ml (1 – 10)	835.TRIO.Y



Code	Description	Volume ul	Division ul	Tip type	Weight gr	% Error - vol .min./medio /max	% Imprecision - vol.min./me dium/max
21825002	Micropipette mod. ACURA 825 from 0.1 - 2 µl	0.1-2	0.002	ultra 10 ul	88	<±6.0/<±4.0/<± 2.0	<5.0/<3.3/<1.5