

Micropipetta ACURA Electro 926 50 - 1000 μ l

Acura Electro is a successful line that makes electronic pipetting easier and safer than ever. The instrument, controlled by a microprocessor, includes intuitive software that allows you to avoid complicated programming steps. No special skills are required to quickly access the six modes of operation and the different volumes available.

Volume range from 0.1 μ l to 10 ml.

- Optimal ergonomic handle low weight
- Easy to learn use
- Large display, adaptable to reading on the right or left side
- Quick replacement battery
- Fast loading, extended working autonomy
- Adjustable ejector *, compatible with most tips
- Counter of pipetting cycles
- Interchangeable volumetric modules
- Two year warranty
- * Socorex patent

The instrument is equipped with a microprocessor-controlled motor and a rechargeable battery.

- The short rod length and tapered end improve controllability
- Easy access to microtubes and microplates
- Optional longer volumetric modules
- User-friendly volume conversion thanks to interchangeable volume modules.
- Unmatched stability and metrological performance

Operating modes:

- Ordinary pipetting
- Reverse pipetting
- Repetitive dispensing
- Dilution
- Tactile pipetting (titration, gel loading)
- Mixing

START PACKAGE

Recommended when purchasing an Acura[®] electro for the first time, each starter pack contains all the components that allow immediate pipetting under optimal conditions.

The additional pipettes can then be purchased on their own and loaded onto an existing stand.

The starter package includes:

- 1x Acura[®] electro 926 XS electronic micropipette
- 1x Qualitix[®] pipette tip samples
- 1x Individual QC Certificate
- 1x Instructions for use
- 1x Charging Stand
- 1x Quick Power Supply with Power Cord and AC Adapter
- 1x Extra battery pack

ACURA ELECTRO ALONE



0.1 - 2 μ L
0.5 - 10 μ L
1 - 20 μ L
2.5 - 50 μ L
5 - 100 μ L
10 - 200 μ L
50 - 1000 μ L



For further purchases of micropipettes, once the START PACKAGE has been purchased, the Acura Electro line can be purchased as a single micropipette with battery back (without charging stand and without battery charger).
The codes are the same as for "Start Packages" without the final "E".



Code	Description	Volume ul	Division ul	Tip type	% Error - vol.m in./mediate/m ax	% Imprecision - vol.min./med ium/max
21926000E	Micropipetta ACURA Electro 926 50 - 1000 µl	50-1000	5.0	1000 ul	<±0,8/<±0,5/<±0, 4	<0,4/<0,15/<0,1