

ACURA Electro 956 12-channel micropipette 20 - 350 μ l

Acura Electro is a successful line that makes electronic pipetting easier and safer than ever.

The instrument, controlled by a microprocessor, includes intuitive software that allows you to avoid complicated programming steps. No special skills are required to quickly access the six modes of operation and the different volumes available.

- Optimal ergonomic handle low weight
 - Easy to learn use
 - Large display, adaptable to reading on the right or left side
 - Quick replacement battery
 - Fast loading, extended working autonomy
 - Adjustable ejector *, compatible with most tips
 - Counter of pipetting cycles
 - Interchangeable volumetric modules
- * Socorex patent

Two year warranty.

The instrument includes a microprocessor controlled motor and a rechargeable battery.

Wide choice of 8 and 12 channel models

- Low weight and perfect grip
- Sequential ejection of the tip
- The body is also suitable for all single-channel volumetric modules

Eight models cover volumes between 0.5 and 350 μ l.

Other volumetric modules available separately.

Operating modes:

- Ordinary pipetting
- Reverse pipetting
- Repetitive dispensing
- Dilution
- Tactile pipetting (titration, gel loading)
- Mixing

START PACKAGE

Recommended when purchasing an Acura[®] electro for the first time, each starter pack contains all the components that allow immediate pipetting under optimal conditions. The additional pipettes can then be purchased on their own and loaded onto an existing stand.

The starter package includes:

- 1x Acura[®] electro 956 electronic micropipette
- 1x Qualitix[®] pipette tip samples
- 1x Individual QC Certificate
- 1x Instructions for use
- 1x Charging Stand
- 1x Quick Power Supply with Cable and AC Adapter
- 1x Extra battery pack



added
NEW
features

0.5 - 10 μ L
2.5 - 50 μ L
10 - 200 μ L
20 - 350 μ L



ACURA ELECTRO ALONE

For further purchases of micropipettes, once the START PACKAGE has been purchased, the Acura Electro line can be purchased as a single micropipette with battery back (without charging stand and without battery charger).

The codes are the same as the "Start Package" without the final "E" suffix.



Code	Description	Volume ul	Division ul	Tip type	% Error - vol.m in./mediate/m - vol.min./medium/max	% Imprecision - vol.min./medium/max
21956352E	ACURA Electro 956 12-channel micropipette 20 - 350 µl	20-350	5.0	350 ul	<±1,0/<±0,8/<±0,6	<0,6/<0,4/<0,25