

## Sodium heparinated microhematocrit capillaries 80 iu-ml vol. ul 75

Microhematocrit glass capillaries are used for centrifugation of blood samples. Heparin capillaries, delaying blood clotting, should be used when samples are not immediately waxed and centrifuged.

- made of high quality glass
- in accordance with DIN ISO 12 772 standards
- length about 75 mm  $\pm$  0.5 mm
- thickness: about 0.2  $\pm$  0.025mm
- sodium heparin (80 iu / ml  $\pm$  30%)
- for single use
- approved for in vitro diagnostic (IVD) applications, according to IVD directive 98/79 / EC, with CE mark.



Packaging: dispenser of 100 pieces, packaging 10 dispensers of 100 pieces

Cnd: W0501010203

Code	Description	internal measures mm diam x length	external dimensions mm diam x length	Volume ul	Pieces per pack	Packing
1159001	Sodium heparinated microhematocrit capillaries 80 iu- ml vol. ul 75	1,1-1,2 x 75	1,5-1,6 x 75	75	100	10x100