

## Swab rod with aluminum support and cotton tip

### sterile in single pack

Sterile cotton / viscose examination swabs on wooden / plastic / aluminum support.  
Outpatient use. Single use. LATEX FREE



#### STORAGE INSTRUCTIONS.

Store in a cool, dry place at room temperature and away from heat sources. Do not expose to direct sunlight, UV light and fluorescent lamps.  
If the packaging is damaged or wet, discard the product. Do not use beyond the date of expiry indicated on the packaging.

#### USEFUL LIFE OF THE PRODUCT

5 years

#### PRECAUTIONS AND SAFETY

Warning: the product is combustible. Combustion can produce toxic fumes: carbon monoxide carbon, carbon monoxide, organic acids. Extinguishing media: water, powder, CO<sub>2</sub>.

#### COMPOSITION

Cotton wool / Viscose  
PVC / birch wood / aluminum

#### STERILIZATION

Through ETO. The sterilization method has been validated in accordance with the UNI EN ISO 11135 standards.

#### METHOD OF DISPOSAL

To be disposed of as hospital waste according to current legislation.

Compliant with: EEC Directive 93/42, Directive 47/2007, UNI EN ISO 11135, EC Regulation 1907/2006 REACH

Class I medical device sterile pursuant to Legislative Decree 37/2010 in implementation of DIR. 93/42 EC smi

| Code     | Description                                   | Rod material | Buffer material | Length mm | Sterile         | Pieces per pack | Packing |
|----------|---|--------------|-----------------|-----------|-----------------|-----------------|---------|
| 05417001 | Swab rod with aluminum support and cotton tip | aluminum     | cotton          | 150       | YES single pack | 500             | 6X500   |