

Container for histology 500 ml - formalin 300 ml

pre-filled with 300 ml of 10% neutral buffered formalin

Universal fixative packaged in hermetically sealed containers screen-printed or labeled, in compliance with current legislation on classification, labeling and packaging of dangerous preparations.



The pre-filled containers represent the optimal solution for the safety of the operator and for the transport of histological samples.

PP container.

Watertight blue PE screw cap for containers from 10ml to 125ml.

Hermetic pressure lid for containers from 250 to 5000 ml.

PVC label resistant to wear and water, alcohol, solvents.

Scratch resistant ink resistant to water and alcohol.

After the first opening, the reagent must be considered valid until the expiry date indicated as long as it is properly stored.

Shelf life of the product: 2 years.

CE - IVD mark.

CND W01030705

Code	Description	Cap.nomin.ml.	Pieces per pack.
05010500	Container for histology 500 ml - formalin 300 ml	300	6