



Over 50 years in cooperation with the health world, the scientific research, the environmental and quality control world

Conductivity calibration solution according to EN 10204

ref. temp. + 25°C - 147μ S/cm - 250 ml

Conductivity calibration solution according to EN 10204.

Conductive solutions are used for the calibration of most conductivity instruments.

To calibrate an instrument, choose a conductive solution closest to the measurement range and calibrate the instrument for optimal system accuracy.

Victor | Vic

Available with different ranges

Available in different formats

Calibration certificate is available on request

Code	Description	conductivity at 25°C	Quantity ml
7506022	Conductivity calibration	147 μS/cm	250
	solution according to EN 10204		

1 / 1