



Over 50 years in cooperation with the health world, the scientific research, the environmental and quality control world

## Conductivity calibration solution according to EN 10204

## ref. temp. + 25°C - $1413 \mu S/cm$ - 500 ml

Conductivity calibration solution according to EN 10204.

Conductive solutions are used for the calibration of most conductivity instruments.

To calibrate an instrument, choose a conductive solution closest to the measurement range and calibrate the instrument for optimal system accuracy.

World Carlo Carlo

Available with different ranges

Available in different formats

Calibration certificate is available on request

Code	Description	conductivity at 25°C	Quantity ml
7506031	Conductivity calibration	1413 μS/cm	500
	solution according to EN 10204		

1/1