



Over 50 years in cooperation with the health world, the scientific research, the environmental and quality control world

Conductivity calibration solution according to EN 10204

ref. temp. + 25° C - 84μ S/cm - 500 mI

Conductivity calibration solution according to EN 10204.

Conductive solutions are used for the calibration of most conductivity instruments.

To calibrate an instrument, choose a conductive solution closest to the measurement range and calibrate the instrument for optimal system accuracy.

Available with different ranges

Available in different formats

Calibration certificate is available on request



Code	Description	conductivity at 25°C	Quantity ml	
7506033	Conductivity calibration	84 μS/cm	500	
	solution according to EN 10204			

1/1