

IoT-Line Precision Balance KERN 572













· Protective working cover included with delivery

Technical data

- · Large backlit LCD display, Digit height 21 mm
- Dimensions weighing surface, stainless steel A Ø 106 mm B Ø 150 mm
- C W×D 160×200 mm, see larger picture
- Overall dimensions W×D×H 180×310×85 mm
- · Net weight A, B approx. 2,4 kg C approx. 2,8 kg
- Permissible ambient temperature -10 °C/40 °C

All-rounder, e.g. as precision balance in the laboratory or in harsh industrial applications, ideal for the diverse possibilities of Industry 4.0 applications

Features

- · Thanks to the many typical laboratory functions, such as, for example, recipe function, percentage determination, combined with the high level of precision, the KERN 572 is a reliable partner for day-to-day work in the laboratory
- · The robust version, typical industrial functions, such as piece-counting, vibrationfree weighing and the large weighing ranges also make these balances ideal for all industrial applications, where a high level of precision is required
- · KERN Universal Port (KUP): permits the connection of an external KUP interface adapter, such as, for example, RS-232, USB, Bluetooth, WiFi or Ethernet, for the exchange of data and control commands, without any installation outlay
- · KERN Communication Protocol (KCP): The KCP permits searching and remote control of the ba-

- lance using external control devices or computers
- · For further information on KUP and KCP see page 20/21
- · Standardised, simplified concept of operation
- · Freely programmable weighing unit, e.g. display direct in special units such as length of thread g/m, paper weight g/m², or similar
- The robust aluminium diecast housing maintains the stability, protects the weighing technology elements and is robust enough to cope with everyday use
- · A Ring-shaped draught shield standard, only for models with weighing plate size A, weighing space Ø×H 157×43 mm
- Level indicator to level the balance precisely
- Underfloor weighing: load support on the underside of the balance for models. Loop for underfloor weighing, standard for models with weighing plate size A

Accessories

- · Protective working cover, scope of delivery: 5 items, KERN 572-A02S05
- · Internal rechargeable battery pack, operating time up to 48 h without backlight, charging time approx. 8 h, KERN YKR-01
- · External data interface RS-232, interface cable included, KERN KUP-01
- · External data interface USB, interface cable included, KERN KUP-03
- · Interface adapter Ethernet, KERN KUP-04
- Interface adapter WiFi, KERN KUP-05
- · Bluetooth interface adapter, KERN KUP-06
- Extension box for connecting up to three interfaces in parallel, KERN KUP-13
- B, C Hook for underfloor weighing, KERN 572-A03
- 1 Large glass draught shield with 3 sliding doors for easy access to the items being weighed. Weighing space W×D×H 150×140×130 mm, for models with weighing plate size A, KERN 572-A05























































	, , , , , , , , , , , , , , , , , , ,						ЪС
Model	Weighing capacity [Max]	Readability [d]	Reproducibility	Linearity	Resolution	Weighing plate	Options DAkkS Calibr. Certificate DAkkS
KERN	g	g	g	g	Points		KERN
572-30	240	0,001	0,001	± 0,003	240.000	Α	963-127
572-31	300	0,001	0,002	± 0,005	300.000	Α	963-127
572-32	420	0,001	0,002	± 0,005	420.000	Α	963-127
572-33	1600	0,01	0,01	± 0,03	160.000	В	963-127
572-35	2400	0,01	0,01	± 0,03	240.000	В	963-127
572-37	3000	0,01	0,02	± 0,05	300.000	В	963-127
572-39	4200	0,01	0,02	± 0,05	420.000	В	963-127
572-45	12000	0,05	0,05	± 0,15	240.000	С	963-128
572-55	20000	0,05	0,1	± 0,25	400.000	С	963-128
572-43	10000	0,1	0,1	± 0,3	100.000	С	963-128
572-49	16000	0,1	0,1	± 0,3	160.000	С	963-128
572-57	24000	0,1	0,1	± 0,3	240.000	С	963-128





Internal adjusting Quick setting up of the balance's accuracy with internal adjusting weight (motordriven)



Adjusting program CAL For quick setting up of the balance's accuracy. External adjusting weight required



EasyTouch Suitable for the connection, data transmission and control through PC or tablet



Memory

Balance memory capacity, e.g. for article data, weights, PLU etc.



Alibi memory Secure, electronic archiving of weighing

results, complying with the 2014/31/EU standard.



KERN Universal Port (KUP)

allows the connection of external KUP interface adapters, e.g. RS-232, RS-485, SB, Bluetooth, WIFI, Analogue, Ethernet etc. for the exchange of data and control commands, without installation effort



RS-232 Data interface

To connect the balance to a printer, PC or network



RS-485 Data interface

To connect the balance to a printer, PC or other peripherals. Suitable for data transfer over large distances. Network in bus topology is possible



USB Data interface

To connect the balance to a printer, PC or other peripherals



Bluetooth* Data interface

To transfer data from the balance to a printer, PC or other peripherals



WIFI Data interface

To transfer data from the balance to a printer, PC or other peripherals



Control outputs

(optocoupler, digital I/O) To connect relays, signal lamps, valves, etc.



Analogue interface to connect a suitable

peripheral device for analogue processing of the measurements



Interface for second balance

For direct connection of a second balance



Network interface

For connecting the scale to an Ethernet network



KERN Communication Protocol (KCP)

It is a standardized interface command set for KERN balances and other instruments, which allows retrieving and controlling all relevant parameters and functions of the device. KERN devices featuring KCP are thus easily integrated with computers, industrial controllers and other digital systems



GLP/ISO log intern

The balance displays weight, date and time, independent of a printer connection



GLP/ISO log Printer

With weight, date and time. Only with KERN printers.



Piece counting

Reference quantities selectable. Display can be switched from piece to weight



Recipe level A

The weights of the recipe ingredients can be added together and the total weight of the recipe can be printed out



Recipe level B

Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through display



Totalising level A

The weights of similar items can be added together and the total can be printed out





Percentage determination Determining the deviation

in % from the target value (100%)



Weighing units

Can be switched to e.g. nonmetric units. See balance model. Please refer to KERN's website for more details



Weighing with tolerance range (Checkweighing)

Upper and lower limiting can be programmed individually, e.g. for sorting and dosing. The process is supported by an audible or visual signal, see the relevant model



Hold function

(Animal weighing program) When the weighing conditions are unstable, a stable weight is calculated

as an average value



Protection against dust and water splashes IPxx

The type of protection is shown in the pictogram



Suspended weighing Load support with hook on the underside of the balance



Battery operation Ready for battery opera-

tion. The battery type is specified for each device



Rechargeable battery pack Rechargeable set



Universal plug-in power supply

with universal input and optional input socket adapters for A) EU, CH, GB B) EU, CH, GB, US C) EU, CH, GB, US, AUS



Plug-in power supply 230V/50Hz in standard version for EU, CH, or AUS version available



Integrated power supply unit

Integrated in balance. 230V/50Hz standard EU. More standards e.g. GB, USA or AUS on request



Weighing principle Strain gauges

Electrical resistor on an elastic deforming body



Weighing principle Tuning fork

A resonating body is electromagnetically excited, causing it to oscillate



Weighing principle **Electromagnetic force** compensation

Coil inside a permanent magnet. For the most accurate weighings



Weighing principle Single cell technology

Advanced version of the force compensation principle with the highest level of precision



Conformity Assessment The time required for conformity assessment is specified in the pictogram



DAkkS calibration possible (DKD) The time required for

DAkkS calibration is shown in days in the pictogram



Factory calibration (ISO)

The time required for Factory calibration is shown in days in the pictogram



Package shipment

The time required for internal shipping preparations is shown in days in the pictogram



Pallet shipment

The time required for

internal shipping preparations is shown in days in the pictogram



