MILLMIX 20



Bead beater

Discover our high-performance ball mill, the ultimate solution for fast and efficient sample disintegration. Ideal for applications in forensic science, pharmacy, metallurgy, and beyond. With a range of available accessories tailored to different needs, this mill is a must-have for precision and versatility in your lab.

| TECHNICAL SPECIFICATIONS | | |
|---------------------------|---|--|
| Model | Millmix 20 | |
| Max. speed | 30 Hz (3600 strokes/min) | |
| Stroke | 32 mm | |
| Max. sample size | ~20 mm | |
| Max. volume | 2 x 50 ml | |
| Max. recomended capacity | up to 2 x 600 g | |
| Display | LED | |
| Timer | 5s to 99 min 50s and continuous, step 1s, 10s | |
| Max. power consumption | 200 W | |
| Dimension (H x W x D) | 240 x 385 x 420 mm | |
| Height (with open lid) | 465 mm | |
| Crating (H x W x D) | 300 x 500 x 500 mm | |
| Weight (Net/Crated) | 42,0 / 58,0 kg | |
| Motor type | induction brushless motor | |
| Permitted working ambient | "temperature 5°C to 40°C, | |
| Tommittod Working ambient | rh 85% non-condensing" | |

| ORDER INFORMATION | |
|-------------------|---|
| 230 V, 50/60 Hz | Millmix 20 |
| Code | 670.000 |
| Order ref. | 289728 |
| 120 V, 50/60 Hz | Millmix 20 120 V |
| Code | 1.670.000 |
| Order ref. | 308251 |
| 230 V, 50/60 Hz | Millmix 20 (Stainless steel Pharmaceutical) |
| Code | 670.800 |
| Order ref. | 314117 |





Jars and adapters

| a | | |
|-----------------|----------|-----|
| Stainless steel | | |
| Jar 1,5 ml | set of 2 | |
| Code | 670.210 | 9 |
| Order ref. | 289934 | |
| Jar 5 ml | set of 2 | |
| Code | 670.220 | |
| Order ref. | 289945 | |
| Jar 10 ml | set of 2 | |
| Code | 670.230 | |
| Order ref. | 289954 | |
| Jar 25 ml | set of 2 | E P |
| Code | 670.240 | |
| Order ref. | 289958 | FL |
| Jar 35 ml | set of 2 | E |
| Code | 670.250 | |
| Order ref. | 289963 | F |
| Jar 50 ml | set of 2 | E |
| Code | 670.260 | |
| Order ref. | 289969 | |

| Hard steel | | |
|------------|----------|--|
| Jar 1,5 ml | set of 2 | |
| Code | 670.310 | 9 |
| Order ref. | 293494 | |
| Jar 5 ml | set of 2 | |
| Code | 670.320 | |
| Order ref. | 293523 | |
| Jar 10 ml | set of 2 | |
| Code | 670.330 | |
| Order ref. | 293526 | |
| Jar 25 ml | set of 2 | R.3 |
| Code | 670.340 | |
| Order ref. | 293529 | FL |
| Jar 35 ml | set of 2 | |
| Code | 670.350 | |
| Order ref. | 293542 | The state of the s |
| Jar 50 ml | set of 2 | (C) |
| Code | 670.360 | |
| Order ref. | 293546 | F |
| | | |

| Special applications | | |
|-------------------------------------|----------|-----|
| Jar 25 ml Tungsten Carbide | set of 2 | 6 |
| Code | 670.440 | |
| Order ref. | 314806 | 6 |
| Jar 35 ml Zircon Oxide | set of 2 | Æ |
| Code | 670.470 | |
| Order ref. | 385739 | 100 |
| Jar 25 ml SS - forensic application | set of 2 | E |
| Code | 670.620 | - |
| Order ref. | 331598 | 6 |
| Jar 50 ml POM | set of 2 | E |
| Code | 670.560 | |
| Order ref. | 310357 | 6 |

| 1 | Adapters POM | | |
|---|---------------------------------------|----------|----|
| | POM Adapter 12 x 0,2 ml | set of 2 | |
| | Code | 670.515 | E. |
| | Order ref. | 290086 | |
| | POM Adapter 12 x 2 ml | set of 2 | |
| | Code | 670.510 | 8 |
| | Order ref. | 290081 | |
| | POM Adapter 6 x 2 ml | set of 2 | 0 |
| | Code | 670.500 | C. |
| | Order ref. | 290040 | |
| | POM Adapter 6 x 2 ml - screw cap tube | set of 2 | 0 |
| | Code | 670.505 | 6 |
| | Order ref. | 290042 | |
| | | | |



Domel

Holders

| Holders | | |
|--|----------|--------|
| H3/MTP POM- holder for 3 x standard or 1 deepwell | set of 2 | |
| Code | 670.450 | |
| Order ref. | 289980 | - I |
| H3/MTP SS -Stainless steel holder for 3 x standard or 1 deepwell | set of 2 | |
| Code | 670.455 | |
| Order ref. | 354457 | |
| H 5/50 SS - Holder for 50 ml conical tubes - 5 tubes in holder | set of 2 | RE |
| Code | 670.530 | |
| Order ref. | 314116 | 16 |
| H 4/50 SS- Holder for 50ml conical tubes - 4 tubes in holder | set of 2 | |
| Code | 670.650 | |
| Order ref. | 366861 | B) II |
| H 5/15 SS - Holder for 15 ml conical tubes - 5 tubes in holder | set of 2 | NE |
| Code | 670.540 | |
| Order ref. | 322538 | 120 |
| H 6/15 SS - Holder for 15 ml conical tubes - 6 tubes in holder | set of 2 | |
| Code | 670.660 | |
| Order ref. | 379173 | |
| H 10/5 SS - Holder for 5ml microtubes - 10 tubes in holder | set of 2 | NE B |
| Code | 670.540 | |
| Order ref. | 322538 | |
| H 6/5 SS - Holder for 5ml microtubes - 6 tubes in holder | set of 2 | |
| Code | 670.640 | () all |
| Order ref. | 382297 | 1971 |
| H 24/2 POM SC - Holder for 1,5/2,0ml microtubes with screw cap | set of 2 | |
| Code | 670.680 | |
| Order ref. | 311345 | |
| | | |



Grinding balls



Stainless steel balls

| dia. 3 mm | Set of 20 |
|------------|-----------|
| Order ref. | 311988 |
| dia. 5 mm | Set of 20 |
| Order ref. | 311600 |
| dia. 7 mm | Set of 20 |
| Order ref. | 311989 |
| dia. 10 mm | Set of 10 |
| Order ref. | 311603 |
| dia. 12 mm | Set of 10 |
| Order ref. | 311605 |
| dia. 20 mm | Set of 5 |
| Order ref. | 323420 |
| dia. 25 mm | Set of 5 |
| Order ref. | 311607 |

Hard steel balls

| dia. 3 mm | Set of 20 |
|------------|-----------|
| Order ref. | 311998 |
| dia. 5 mm | Set of 20 |
| Order ref. | 311610 |
| dia. 7 mm | Set of 20 |
| Order ref. | 311999 |
| dia. 10 mm | Set of 10 |
| Order ref. | 311621 |
| dia. 12 mm | Set of 10 |
| Order ref. | 311622 |
| dia. 25 mm | Set of 5 |
| Order ref. | 311623 |

Others

| • ******* | |
|------------------|-----------|
| Agate | |
| dia. 10 mm | Set of 10 |
| Order ref. | 311640 |
| Zircony Oxide | |
| dia. 15 mm | Set of 5 |
| Order ref. | 371033 |
| dia. 20 mm | Set of 5 |
| Order ref. | 371034 |
| Tungsten Carbide | |
| dia 3mm | Set of 20 |
| Order ref. | 315284 |
| dia 5 mm | Set of 20 |
| Order ref. | 315286 |
| dia 7 mm | Set of 20 |
| Order ref. | 315288 |
| dia. 15 mm | Set of 5 |
| Order ref. | 375920 |
| dia. 20 mm | Set of 5 |
| Order ref. | 315289 |
| | |

