

Flask 300 cm², click-vent cap

Absence of heavy metals, cytotoxic agents, DNA / RNA, DNase / RNase, ATP, BSE / TSE.

Pyrogen levels <0.06 EU / ml.

Cap with 0.22 micron PTFE hydrophobic filter. The filter avoids contamination and is water-repellent towards liquids allowing an excellent gas exchange.

Highly uniform activated surface

Extreme transparency of plastic.

Neck geometry that reduces the risk of soil spillage and at the same time guarantees total access to pipettes, scrapers / cell lifters, absence of "dead corners" with consequent recovery of 100% of the cells.

Stable shape, with the possibility of easily stacking multiple flasks.

Bilateral area of clearly visible marking (white rectangle placed on the 2 upper sides of the flask).

Volume graduations shown in white for more precise optical control of fill volume



Click-vent cap : with double position (complete closure or controlled ventilation opening)

Filter cap : equipped with a hydrophobic membrane in PTFE with porosity 0.22 µm

For flasks with a "Vent " type cap:

- A " **click** " and the cap with the word «VENT» is in the «VENT» position for ventilation.

For a visual check : the rectangle placed on the cap, slightly raised, is visible at " **12:00** " . - TURNING the "VENT" **screw cap 90 °** to close the tube and no gas passage / aeration.

For a visual check : the rectangle placed on the Cap, slightly raised, is at " **3:00** ". - See more at: <http://www.delchimica.com/it-IT/articoli/ricerca-61.aspx#sthash.ebuv2wTu.dpuf>

Code	Description	Growth area cm ²	Cap version	Volume ml	Dimensions mm	Pieces per pack.
2790300	Flask 300 cm ² , click-vent cap	300	VENT	1360	275x170x50	18