


Cell scraper 20 x 300 mm

With a slight downward pressure on the handle of the scraper it is possible to vary the angle of the mobile spatula of the scraper itself. The slight rotation of the handle then allows the spatula to rotate in the desired direction.

Mobile spatula:

Free rotation of the tip of the scraper's movable spatula
Full access to all corners
Different sizes of spatula
Small elevations on the handle of the scrapers to facilitate the grip

The scrapers are individually packaged in easy-to-open bags.



Code	Description	Blade dimensions mm	Length mm	Pieces per pack.
2799003	Cell scraper 20 x 300 mm	20	300	150