

Automatic pipette with repetitive dosage Brand mod. HandyStep® S

The HandyStep® S automatic repetitive dosing pipette is ideal for routine work and research in the fields of diagnostics, molecular biology, environmental analysis and many more.

HandyStep® S and PD-Tips work on the principle of positive displacement. Even fluids with high viscosity, density or vapor pressure can therefore be metered with maximum precision. The positive displacement allows you to work without contamination because no aerosols are formed.

- greater resistance to chemicals thanks to the use of innovative plastic materials
- easy insertion of the tip: the PD-Tip is now inserted straight from the bottom
- Volumes between 2 µl and 5 ml
- up to 49 repetitions of the dosage
- only 108 g in weight
- CE IVD compliance

Original BRAND PD TIPS microdispenser syringes available on request



Code	Description
53705100	Automatic pipette with repetitive dosage Brand mod. HandyStep® S