

Drug test - suitable in occupational medicine for economic screening of workers

By simply scanning one of the tests in the device (code 71024576, 71024577), it is possible to obtain a complete analysis on your PC.

- automatic analysis of visual reading tests
- internal quality self-control system
- image acquisition, data management / export, results archive

Main benefits:

- automatic - eliminates variability and subjectivity in the interpretation of results
- precise - guarantee of the result with respect to visual tests
- fast - the scanning process takes 15-30 seconds
- flexible - different formats / menus, light weight, semi-portable, easy installation, compatible with Windows

Usable tests:

71024576 DRUG TEST IN GLASS - 7 drugs + antiadulterants

71024577 DRUG TEST IN GLASS - 8 drugs + antiadulterants



Code	Description	Pieces per pack.
71024570	Drug Reader - device for drug test analysis with software in Italian	1
71024576	Drug test per shot glass for 7 drugs + antiadulterants	25
71024577	Drug test for 8 drugs + antiadulterants	25