

SPECI-SPONGE sterile sampling bags 115x230 ml 532

with sponge for environmental sampling of surfaces

The Whirl-Pak® Speci-Sponge® pouch is designed for surface environmental sampling of work areas, equipment, animal carcasses and any other locations where testing for *Listeria*, *Salmonella*, *E. coli*, coliforms and other pathogens of food origin or anthrax is necessary.

The special dry sponge is free of bactericides and has been tested not to be inhibitory. It measures approximately 3.8 x 7.6 x 1.5cm thick when wet, and is sterilized inside a Whirl-Pak® pouch.

It is recommended to wear a sterile glove when cleaning the sponge to prevent contamination.



After the sample has been collected and the sponge inserted into the bag, it can be sent to a laboratory for testing.

| Code | Description | Dim. Mm | Cap. Ml | Th. Mm | Writing band | Pieces per pack. |
|----------|--|---------|---------|--------|--------------|------------------|
| 73001245 | SPECI-SPONGE sterile sampling bags 115x230 ml 532 | 115x23 | 532 | 0.064 | SI | 100 |