

Aluminum rocker mod. B175 for rotor mod. S4-175

Aluminum rocker mod. B175 for rotor mod. S4-175


Max capacity 175 ml

To be completed with the specific adapters as shown in the table below

Adapters for the scales mod. B175 (rotor mod. S4-175)			
Code	Description	Total tubes per rotor	Notes for order
B4100602	12-place adapter mod. I12-2 for microtubes type Eppendorf 1,5 / 2,0 ml (diam. mm 13 x L 75 and 13 x L 100)	48 test tubes	4-piece set
B4100612	7-seat adapter mod. I7-5 / 7 for tubes type VAC 5-7 ml (diam. mm 13 x L 75 and diam. 13 x L 100)	28 test tubes	4-piece set
B4100622	7-seat adapter mod. I7-10 for test tubes type VAC 10 ml (diam. mm 16 x L 100)	28 test tubes	4-piece set
B4100632	4-seat adapter mod. I4-15F for tubes with conical bottom type Falcon 15 ml (diam. mm 17 x L 120)	16 test tubes	4-piece set
B4100642	4-seat adapter mod. I4-15R for round bottom tubes 15 ml	16 test tubes	4-piece set




(diam. mm 17 x L 110)



B4100652	1 seat adapter mod. I1-50F for Falcon type conical bottom tubes 50 ml (diam.mm 29 x L 116)	4 test tubes	4-piece set
B4100662	1 seat adapter mod. I1-50R for round bottom tubes 50 ml (diam.mm 30 x L 110)	4 test tubes	4-piece set
B4100672	1 seat adapter mod. I1-100 for round bottom tubes 100 ml (diam.mm 46 x L 102)	4 test tubes	4-piece set

Code	Description
B4101512	Aluminum rocker mod. B175 for rotor mod. S4-175

Accessories

Foto	Code	Description
	B4101522	