

## Tubes with Na CITRATE for ESR

### 2 ml - Pink cap - flat bottom - 16x60 mm

Snap cap hematology tubes with sodium citrate liquid anticoagulant for ESR.

Tubes with 0.4 ml of sodium citrate (3.8%) for 1.6 ml of blood (VES)

Test tubes available in different sizes Suitable for analysis in hematology and immunohematology.

With fill level indicator label.

Perfect traceability.

Unbreakable.

Particularly transparent medical Polypropylene (PP) test tube.

Pressure cap in colored polyethylene perfectly sealed.

Latex free device.

IVD CE certificate



Code	Description	Nominal volume mL	Diam.x Height mm	Bottom	Cap Color	Conf.	Sub Conf.
05055014	Tubes with Na CITRATE for ESR	0,4 ml of sodium for 1,6 ml of blood	16x60	Flat	Pink	1500	30x50