

Porcelain mortars, with raw inside and enamelled outside, with spout

5 ml

Tissue grinder is composed by two components to be ordered separately: Borosilicate glass mortar (items 6102 to 6150) and PTFE pestle.

Ideal for small amounts of tissue.

PTFE pestle has 0.15 - 0.25 mm clearance with borosilicate mortar.

Volumes stated are working volumes.

Avoid excessive generated heat as could cause binding.

Material: BOROSILICATE GLASS



Code	Description	Internal diameter mm	Height mm	Capacity ml	Pieces per pack.
03106105	Porcelain mortars, with raw inside and enamelled outside, with spout	5	12	135	1