

CAMERA CONTAGLOBULI MALASSEZ

Le camere contaglobuli si differenziano per le griglie di conteggio e per la profondità delle camere. Linee perpendicolari, incise nella superficie della base della camera, formano una griglia che diventa visibile al microscopio.

Approvate per applicazioni (IVD) diagnostici in vitro, secondo la direttiva IVD 98/79 / CE, con marchio CE

A richiesta sono disponibili camere contaglobuli con reticolo "Bright Lined" che facilitano la valutazione del campione grazie alle linee luminose che contrastano bene con lo sfondo scuro.


La profondità della camera Malassez è di 0,2 mm.

La griglia di conteggio copre 2 x 2,5 mm².

I rettangoli hanno un'area di 0,25 x 0,20 = 0,05 mm².

Ciascuno di essi è suddiviso in 20 piccoli quadrati con una superficie di 0,05 mm x 0,05 mm = 0,0025 mm² ciascuno.


Questa camera viene utilizzata per il conteggio delle cellule nel liquido cerebrospinale e per il conteggio dei nematodi.



Malassez

The depth of the Malassez counting chamber is 0.2 mm. The counting grid covers 2 x 2.5 mm². The large rectangles have an area of 0.25 x 0.20 = 0.05 mm². Each of them is subdivided into 20 small squares with an area of each 0.05 mm x 0.05 mm = 0.0025 mm².

This counting chamber is used for counting cells in liquor (cerebro-spinal fluid) and for counting nematodes for example.



Codice	Descrizione	Profondità della camera mm
1221016	Camera conteggio Malassez doppia senza pinze	0,2
1221018	Camera conteggio Malassez doppia con pinze	0,2