


## *Tissue Culture Dish 40 – 150 cm<sup>2</sup>*



Tissue culture dishes from TPP stand for large growth area and secure handling.

### Features of the Tissue Culture Dish:


- TPP surface activation of the growth area optimally enhances the proliferation of the cells.
- The gripping ring: often copied, but never like the original!
- Serrated edges in the gripping ring prevent accidental dropping of the dish bottom
- Lateral yellow inscription field on the lid
- Stacking ring on the lid and corresponding on the base result in an extremely secure stacking feature

- Clock numbering system in the base of the dish divides it into quadrants for quick orientation
- Spacer on the inside of the dish lid for a constant movement of air
- Vents in the dish base for airflow between stacked dishes for consistent growth conditions and prevention of condensation

Impeccable optical clarity enables a visual control of the cells in stacked dishes using transmitted light. This is the basis for the documentation with microphotography.

50 anni in collaborazione con il mondo della sanita', della ricerca scientifica, del controllo ambientale e di qualita'

- 1 Convenient peel-off bags.
- 2 Clock-numbering system (3, 6, 9 and 12) inlaid into the base of the dish allows identification, region of interest definition and documentation of the zones.



Tissue culture dishes are designed for manual handling only.

- 1 Stacked tissue culture dishes.
- 2 The gripping ring: often copied, but never like the original!

The yellow inscription field on the lid and the matte inscription field on the base enable a defined positioning of the lid. The side walls of the dish are not treated for tissue growth.

| Type | Product-No. | Growth surface<br>cm <sup>2</sup> | Internal-<br>Ø<br>mm | Dimensions<br>h<br>mm | Material | Qty /<br>Bag | Qty /<br>Case |
|------|-------------|-----------------------------------|----------------------|-----------------------|----------|--------------|---------------|
| ○    | 93040       | 9.2                               | 34                   | 40 x 11               | PS       | 20           | 900           |
| ○    | 93060       | 22.1                              | 53                   | 60 x 16               | PS       | 14           | 840           |
| ○    | 93100       | 60.1                              | 87                   | 96 x 21               | PS       | 10           | 240           |
| ○    | 93150       | 147.8                             | 137                  | 146 x 21              | PS       | 5            | 100           |



- 3 Lid with yellow inscription field as well as spacers on the inside of the lid.